
R2CP Interface to Sedecal Room Platform (abridged version) 1.4

4. **R2CP**

- 1. Based on ethernet sockets
 - Smart Hub
 - Message format
 - Routing & Filtering
- 2. R2CP SW interface with the room: a new paradigm
 - Workflow oriented vs device oriented
 - Procedure definition
 - Procedure activation
 - System messages

1. Ethernet

- Smart Hub (SH)
 - SW module that runs on Windows
 - Minimizes socket connections on Positioner and Generator
 - All nodes connected to SH
 - SH routes messages between nodes



Message format

Priority	SEQ	Destination Node Number	Issuing Node Number	Index	Subindex	Function	Data Length
2	6	8	8	8	8	8	16


PP	Priority		
SEQ	Sequence Number, local to each issuing number		
DNN	Destination Node Number O Broadcast mesage (ANSWER/EVENT) NN Message to node number NN (SET or GET)		
INN	Issuing Node Number		
Index	 Parameter Group 0x20 Generator 0x40 Positioner 0x60 Dosimeter 0x70 Custom 0x80 GUI & Fileshare (System Messages) 0xA0 Network 0xE0 Patient Workflow 0xF0 Message Workflow 		
Subindex	Parameter		

Function	SETGETANSWER/EVENTCOMMAND PROCESSED
Data Length	Number of bytes in the data field


Routing and filtering: rules

Implicit Rule	All messages received from any node with DNN = NN will be forwarded to node number NN	
Broadcast Message	Smart Hub will forward to NN node all the messages with DNN=0x00 when NN node is registered to receive messages from the INN and the range of Index/Subindex specified in the rules	
	Rule 1 INN=A, specifies the issuing node number the rule applies to. Index Subindexes: they can be defined as individual indexes and ranges.	
	Rule 2 INN=B Index Subindexes	

- Routing & Filtering: Nested SH's: If destination node is not in the list of local nodes, message will be forwarded to the SH in which node is local.
 - Two positioners & 2-tube generator.



 Adding an additional console is as simple as connecting two instances of SH with the same local nodes connected (different NN).



- R2CP SW interface with the room: a new paradigm
 - Workflow oriented vs Device oriented
 - Related settings defined at once
 - Unified reporting of system events and messages to the user
 - Native support for RF

Intuitive

"Perform a table top procedure for abdomen. Use a wifi detector"



Procedure

- Standard RAD
- 1 Exposure: 1 set of exposure parameters
- 1 Position: table.
- 1 Filter
- 1 Collimator FOV
- Workstation → Place → Table top → Image receptor over the table
- Image receptor: 43x43 wifi detector

Implementing worflow:

- Procedure definition
- Load DB's
- Assign DB's
- Activate DB's

Define Load Assign Activate

- Workstation
 - Refers to image receptor location and the devices involved for the xray procedure to be performed

DATA	FORMAT		
Workstation ID	131		
Description	String		
Generator #	Generator Node Number		
Tube #	Tube number		
Positioner #	Positioner Node Number		
Image Receptor ID	1255		
Image Receptor	0 – Free		
Location	1 – Table		
	2 – Wallstand		
AEC Chamber	0 – No AEC		
Connected to	1 – Generator		
	2 – eBox		
AEC Input	0 - No AEC		
	14 – AEC Input		

Procedure

- Type of XRAY examination: generator procedure type and positioner procedure type.
- Workstation
- Number of exposures
- Number of positions filters collimator FOV's

BYTE	DATA	FORMAT
1	Procedure ID	[1255]
2	Generator Procedure Type	0 – Not defined
		1 – Std RAD
		2 – Stitching
		3 – Tomography
		4 – DSI: single energy, multi energy, tomosynthesis.
		6 – CINE
		7 – DSA
		100 – Std Fluoro
		101 – Roadmap Mask Fluoro
		102 – Roadmap Fluoro
3	Positioner Procedure Type	0 – Not defined
		1 – Std Position
		2 – Stitching
		3 – Tomography
4	Handswitch / Footswitch ID	4 - Tomosynthesis [1255]
4	Handswitch / Footswitch id	[1233]
5	Activate When	0, Do not Activate procedure when Footswitch/Handswitch pressed. Procedure is
	Handswitch/Footswitch is	activated with Patient Workflow message.
	pressed	
		1, Activate procedure when Footswitch/Handswitch pressed
6	Workstation	[1255]

7	Total Number of Exposure Data Banks in Procedure	[1255]		
8	Total Number of Positioning	[1255]		
	Data Banks in Procedure			
9	Total Number of Collimator FOV	[1255]		
	Data Banks in Procedure			
10	Total Number of Collimator Filter	[1255]		
	Data Banks in Procedure			
11	Generator Data Bank	0 – N/A		
	Sequencing	1 – Activate next Exposure Data Bank		
		2 - Pause Procedure until next Exposure Data Bank is activated by a Patient		
		Workflow message		
12	Positioning Data Bank	0 – N/A		
	Sequencing	1 - Activate next Positioning Data Bank		
		2 - Pause Procedure until next Positioning Data Bank is activated by a Patient		
		Workflow message		
13	Collimator FOV Data Bank	0 – N/A		
	Sequencing	1 – Activate next Collimator FOV Data Bank		
		2 - Pause Procedure until next Collimator FOV Data Bank is activated by a Patient		
		Workflow message		
174	Collimator Filter Data bank	0 – N/A		
	Sequencing	1 – Activate next Collimator Filter Data Bank		
		2 – Pause Procedure until next Collimator Filter Data Bank is activated by a Patient Workflow message		

Data Banks (DB's): store the values required by the procedure

- RAD or FL: contains RAD or FL exposure parameters
- Position: contains position information
 - Parameters that define a position
 - Index to a position from library
- Filter: contains filter information
- Collimator FOV: trans and long collimator blades opening

Procedure definition: R2CP Patient Workflow - SET Procedure

 Message sent from image acquisition workstation to generator and positioner

Image Acquisition Workstation

Procedure ID = 1

Generator Procedure Type = Std RAD

Positioner Procedure Type = Std Position

Number of Exposure DB's = 1

Number of Positioning DB's = 1

Number of Filter DB's = 1

Number of Collimator FOV DB's = 1

Generator

Procedure ID = 1

Generator Procedure Type = Std RAD

Positioner Procedure Type = Std Position

Number of Exposure DB's = 1

Positioner

Procedure ID = 1

Generator Procedure Type = Std RAD

Positioner Procedure Type = Std Position

Number of Positioning DB's = 1

Number of Filter DB's = 1

Number of Collimator FOV DB's = 1

Load RAD exposure DB: R2CP Generator message from image acquisition workstation to generator

BYTE	DATA	FORMAT
1	Data Bank ID	1
2(1-4)	Patient Size	0: Small
		1: Medium
		2: Large
2(5-8)	Pediatric	0: Non Pediatric
		1: Pediatric
3(1-6)	Technique Mode	0: 0Point 1: 1Point 2: 2Point 3: 3Point
4 5	kV * 10	[20150] * 10
6 7 8	mAs * 1000	[0.13200] * 1000

9	mA * 100	[11000] * 100
10		
11		
12	ms * 100	[120000] * 100
13		
14		
15	Minimum	[120000]
16	Integration	
	Time (ms)	
17	Maximum	[120000]
18	Integration	
	Time (ms)	
19	Focal spot	0: small
		1: large
		2: super small
20	AEC	0: low
	Sensitivity	1: medium
		2: high
21	AEC Density	[-8 8] (INT8)

22	AEC Chamber off(0) / on(1) status	Bit 1: left AEC status Bit 2: center AEC status Bit 3: right AEC status Bit 4: landscape (0) / portrait (1) AEC orientation
23	Tube Power Limit (%)	[0100]
24 25	FPS x 10	[01200] 0 – Single Shot
26	Tracking ID	[0255] 0 – None
27(1)	kV Scan	0 – Not Active 1 – Active
27(2)	mA Scan	0 – Not Active 1 – Active

Load Fluoro exposure DB: R2CP Generator message from image acquisition workstation to generator

ВҮТЕ	DATA	FORMAT
1	Data Bank ID	[1255]
2(1-4)	Patient Size	0: Small
		1: Medium
		2: Large
2(5-8)	Pediatric	0: Non
		Pediatric
		1: Pediatric
3 4	kV * 10	[40125] * 10
5	mA * 100	[0.0130] *
6		100
7 8	ms * 100	[11000] * 100
9		

10	Maximum	[11000]
11	Integration	
	Time (ms)	
12	PPSx10	[01200] 0 – Continuous
12	FFOXIU	[01200] 0 = Continuous
13		
14	ABC	0 ABC OFF
		1 ABC ON
15	High Dose	0 High Dose OFF
		1 High Dose ON
16	KV SCAN	0 KV SCAN OFF
		1 KV SCAN ON
17	QbyPPS	0 OFF
		1 ON
18	Dose Level ID	[0255]
		0 – None
19	Curve ID	[0255]
		0 – None

Load position DB: R2CP Positioner message from image acquisition workstation to positioner

ВҮТЕ	DATA	FORMAT
1	Data Bank ID	1
2	Programmable	18
3	Position ID	

Load filter DB: R2CP Positioner message from image acquisition workstation to positioner

вуте	DATA	FORMAT
1	Data Bank ID	1
2	Filter	2

Load Collimator FOV DB: R2CP Positioner message from image acquisition workstation to positioner

BYTE	DATA	FORMAT	
1	Data Bank ID	1	
2	Collimator Trans 1	200	
4 5	Collimator Trans 2	200	
6 7	Collimator Long 1	200	
8 9	Collimator Long 2	200	
10 11	Max Collimator Trans 1	2150	
12 13	Max Collimator Trans 2	2150	
14 15	Max Collimator Long 1	2150	
16 17	Max Collimator Long 2	2150	

Assign exposure DB 1 to Procedure 1: R2CP Generator message from image acquisition workstation to generator

ВҮТЕ	DATA	FORMAT
1	Procedure ID	1
2	Exposure	1
	Sequence	
	Number	
3	Exposure Data	1
	Bank ID	

Assign position DB 1 to Procedure 1: R2CP Positioner message from image acquisition workstation to positioner

BYTE	DATA	FORMAT	
1	Procedure ID	1	
2	Positioning Program Sequence Number	1	
3	Positioning Data Bank ID	1	

Assign collimator filter DB 1 to Procedure 1: R2CP Positioner message from image acquisition workstation to positioner

BYTE	DATA	FORMAT	
t	Procedure ID	1	
2	Collimator Filter Sequence Number	1	
3	Collimator Filter Data Bank ID	1	

Assign collimator FOV DB 1 to Procedure 1: R2CP Positioner message from image acquisition workstation to positioner

BYTE	DATA	FORMAT
1	Procedure ID	1
2	Collimator FOV Sequence Number	1
3	Collimator FOV Data Bank ID	1

Activate Procedure and DB's: R2CP Patient Workflow message from image acquisition workstation to generator and positioner

вуте	DATA	FORMAT
1	Procedure ID	1
2	Command	Activate Procedure
3	Exposure Data Bank Sequence Number	1 Activate Exposure Data Bank Sequence Number 1
4	Positioning Data Bank Sequence Number	Activate Positioning Data Bank Sequence Number 1
5	Collimator FOV Data Bank Sequence Number	Activate Collimator FOV Data Bank Sequence Number 1
6	Collimator Filter Data Bank Sequence Number	Activate Collimator Filter Data Bank Sequence Number

Procedure and DB's activation

- Generator: sets exposure parameters with DB 1 contents
- Positioner
 - Position DB: sets as target position programmable position 18
 - Collimator Filter DB: sets filter 2, collimator will turn filter disc until filter 2 is selected.
 - Collimator FOV DB: sets 40x40 FOV, collimator will move blades to achieve this FOV.

Application example: Dual Energy

- Two RAD exposure DB's
- 1 Position DB
- 1 Collimator Filter DB
- 1 Collimator FOV DB

Dual Energy procedure definition



Single Shot: Procedure ID = 1

Generator Procedure Type = 1 Std RAD

Positioner Procedure Type =1 - Std
Position

Generator Data Bank


Sequencing = 1 – Activate next Exposure Data Bank

Workstation = 1

#Exposure DB's = 2

Exposure 1 DB

Exposure 2 DB



RAD Exposure DB ID = 2

Positioning DB ID = 1

Collimator FOV DB ID = 1

Collimator Filter DB ID = 1

POSITIONER

Single Shot: Procedure ID = 1

Generator Procedure Type = 4 Single energy, multi energy.

Positioner Procedure Type =1 – Std
Position

Positioning Data Bank Sequencing = 0 N/A

Collimator FOV Data Bank Sequencing = 0 N / A

Collimator Filter Data bank Sequencing = 0 N/A

#Positioning DB = 1

#Collimator FOV DB = 1

#Collimator filter = 1

Positioning 1 DB

Collimator FOV 1 DB

Collimator Filter 1 DB

Application example: 3-image stitching

- 3 RAD exposure DB's
- 3 Position DB
- 1 Collimator Filter DB
- 3 Collimator FOV DB

3-Image stitching procedure definition

GENERATOR

Single Shot: Procedure ID = 1

Generator Procedure Type = 1 Std RAD

Positioner Procedure Type =1 - Std
Position

Generator Data Bank

Sequencing = 1 – Activate next Exposure Data Bank

Workstation = 1

#Exposure DB's = 3

Exposure 1 DB

Exposure 2 DB

Exposure 3 DB

RAD Exposure DB ID = 1

RAD Exposure DB ID = 2

RAD Exposure DB ID = 2

Positioning DB ID = 1

Positioning DB ID = 2

Positioning DB ID = 3

Collimator FOV DB ID = 1

Collimator FOV DB ID = 2

Collimator FOV DB ID = 3

Collimator Filter DB ID = 1

POSITIONER

Single Shot: Procedure ID = 1

Generator Procedure Type = 4 Single energy, multi energy.

Positioner Procedure Type =1 – Std Position

Positioning Data Bank
Sequencing = 1, Activate next
Positioning Data Bank

Collimator FOV Data Bank

Sequencing = 0 N / A

Collimator Filter Data bank
Sequencing = 0 N/A

#Positioning DB = 3

#Collimator FOV DB = 3

#Collimator filter =1

Positioning 1 DB

Positioning 2 DB

Positioning 3 DB

Collimator FOV 1 DB

Collimator FOV 2 DB

Collimator FOV 3 DB

Collimator Filter 1 DB

System Messages

- Unified reporting of system events and messages to the user
 - Door open: message to the user and exposure inhibit.
 Message and inhibit cleared when door closes.
 - AEC error: message to the user. Generator is blocked until error (not enough dose) is acknowledged by the user.
 - Incorrect grid inserted: exposure inhibit and message displayed until correct grid is inserted.
 - Fatal error on filament board: Call service.
- Included in R2CP 0x80 Index
 - Unique System Message ID
 - Custom system messages can be added

System Message Parameters: Category

- Warning: message to the user
- Abort Exposure: if it occurs during and exposure, it can abort it.
- Inhibit RX
 - Generator Phase: define if inhibits exposure but allows prepping to ready.
- Inhibit Movements
- User Action required: insert correct grid

System Message Parameters: Clear Conditions

- Clear after display
- Clear after user confirmation
- Clear by source device
- Cannot be cleared

System Message Parameters: GUI Indications

- Display modal popup window
- Icon



Warning



Inhibit



Error



User action required



Emergency button



Information

R/F Systems

R/F or Dynamic Systems: SW Interface

- More than one procedure defined, typically a fluoro procedure and at least 1 RAD procedure.
- Procedures can be activated by a message of by a handswitch / footswitch being pressed.
- Procedure definition message includes handswitch / footswitch ID and activation type.

ВҮТЕ	DATA	FORMAT		
4	Handswitch / Footswitch ID	[1255]		
	Activate When Handswitch/Footswitch is pressed	0, Do not Activate procedure when Footswitch/Handswitch pressed. Procedure is activated with Patient Workflow message.		
		1, Activate procedure when Footswitch/Handswitch pressed		

- A fluoro procedure is usually defined as default, meaning that it gets activated when no handswitch / footswitch is pressed.
- Only one procedure can be active at a time, collimator filter and FOV can be different between procedures.

Application example: 1 fluoro procedure, 1 single shot procedure, 1 multi RAD procedure. Procedure ID=3 is default procedure

Procedure ID	1	2	3
Generator Procedure Type	4 – DSI: single energy, multi energy,	4 – DSI: single energy, multi energy,	100 – Std fluoro
Positioner Procedure Type	1 – Std Position (FREE)	1 – Std Position (FREE)	1 – Std Position (FREE)
Handswitch / Footswitch ID	1	2	3
Activate When	1	1	1
Handswitch/Footswitch is pressed			
Workstation	1	1	1
Total Number of Exposure Data	1	1	1
Banks in Procedure			
Total Number of Positioning Data	1	1	1
Banks in Procedure			
Total Number of Collimator FOV Data	1	1	1
Banks in Procedure			
Total Number of Collimator Filter	1	1	1
Data Banks in Procedure			
Exposure DB ID	1	2	3
Positioning DB ID	1	1	1
Collimator filter DB ID	1	1	2
Collimator FOV DB ID	1	1	2

GENERATOR

Single Shot: Procedure ID = 1

Generator Procedure Type = 4 DSI

Positioner Procedure Type =1 - Std Position (FREE)

Handswitch/Footswitch = 1

Activate when handswitch/footswitch pressed = 1

Workstation = 1

#Exposure DB's = 1

Exposure 1 DB

POSITIONER

Single Shot: Procedure ID = 1

Generator Procedure Type = 4 Single energy, multi energy.

Positioner Procedure Type =1 - Std Position (FREE)

Handswitch/Footswitch = 1

Activate when handswitch/footswitch pressed = 1

Workstation = 1

#Positioning DB = 1

#Collimator filter DB = 1

#Collimator filter FOV = 1

Positioning 1 DB

Collimator Filter 1 DB

Collimator FOV 1 DB

RAD Exposure DB ID = 1

kVp

mA

ms mAs

FPS = 0

Positioning DB ID = 1

Collimator Filter DB ID = 1

Collimator FOV DB ID = 1